Corex™ II REF 300-090 (0100860413000109) REF 300-09C (0100860413000147) REF 300-070 (0100860413000116) REF 300-07C (0100860413000154)

(English) Instructions For Use

I. INTENDED USE

The Corex™ Bone Harvester (Corex) is a supplied sterile, single patient use, manually operated trephine intended for harvesting cancellous bone from various skeletal sites. Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.


Warning: Device is a single patient use item reuse may result in failure of proper actuation and or biologic contamination. These consequences could result in adverse patient effects. This device has sharp cutting surfaces and potentially represents a hazard to the user, patient, and disposer. All appropriate cautious handling should be exercised.

Relative Contraindication: Osteopenic/Osteoporotic patients may have insufficient density of bone to facilitate use of the Corex™ Bone Harvester.

Warning: This device has not been evaluated for infants.

II. CORTICAL DEFECT PREPARATION

The Corex™ is supplied with a removable trocar tip, for the purpose of creating a cortical aperture/bore. If the surgeon user opts to use a cannulated trocar version of the Corex™ Bone Harvester, the cannulated trocar tip should be utilized over a guide pin of approximately 2mm in diameter, which has been advanced through the core of the bone into which the harvester is intended to be utilized. The tip is associated with the distal shaft of the harvester. Carefully remove the protective elastomeric cap, then apply the sharp, trocar tip to cortical bone surface, beneath which one intends to remove cancellous bone.

While applying sufficient but not excessive pressure against the cortical bone, rotate the Corex device by gripping the proximal handle and rotating clockwise and then counter clockwise beyond 45 degrees in either direction. Adjust pressure as needed to control depth of penetration. Avoid plunging trocar beyond proximal collared surface (immediately proximal to sharp trocar edges). Once cortical hole is developed, remove the trocar tip from the Corex, by replacing the elastomeric cap, and either depressing the tab engagement while pulling it off, or pushing the trocar tip off with the supplied removal tool.

In the event that the trocar tip is not used, or is insufficient in developing a cortical defect, use an Osteotome, Capner gauge, drill, etc. Make a cortical window large enough to allow the Corex bone harvester to be redirected along various axes upon successive passes.


III. PRECAUTIONS

Once the trocar tip is removed, Corex's trephine (hollow tube section) is not intended to further penetrate cortical bone. There are circumstances of use when cortical penetration might occur, or is even likely to occur. The surgeon must be keenly aware of this potential hazard, as serious or life threatening injuries to nerves, vessels, or other soft tissue structures external to the bone's cortical margins may result. The surgeon should understand the orientation and depth of penetration associated with the Corex cutting tip, relative to the skeletal anatomy within which it is being applied. Every effort to ensure that the cutting tip remains within the cancellous portion of the bone should be exercised. If there is any doubt, the device's depth of entry and axial orientation should be measured and noted and, if appropriate, the device should immediately be withdrawn from the bone entry site.

Warning: Extreme caution should be utilized when placing the sharp tip of the device near vulnerable, "at risk" structures.


IV. PRIOR TO HARVEST

- 1. A plastic extrusion ramp is provided in the bubble pack. DO NOT DISCARD TAMP- tamp will be used to extrude cancellous bone from device.

V. TO HARVEST CANCELLOUS BONE

- 1. Prior to harvest ensure that the upper handle is rotated fully counter-clockwise to the lower handle.
2. An audible click will indicate harvester is locked in the open position and ready for harvesting.

- 3. Rotate Corex™ in an oscillating clockwise-counter-clockwise manner while simultaneously applying advancing pressure to the harvest bed, by pushing on the proximal handle.
4. When full depth of harvest has been achieved, fully lock device by turning upper handle clockwise relative to the stabilized lower handle (locked position).
5. Rotate locked Corex in a clock-wise manner while gradually withdrawing from harvest site.
6. After withdrawing from harvest site:
- Open Corex completely by turning upper handle in counter-clockwise direction relative to the stabilized lower handle to the open locked Corex position (see "Prior To Harvest" instruction #1).
- Extrude cancellous bone from Corex into sterilized dish using plastic tamp.
- Extrude after each pass to prevent cylindrical graft material from binding inside retaining harvesting cylinder.

SUGGESTIONS FOR ILLIAC CREST HARVESTING

POSTERIOR APPROACH

[SUGGESTED SIZE- 9MM] 300-090 and 300-09C

- Harvest from large dilated area of PSIS (Poster Superior Iliac Spinos Process).

Warning: BE MINDFUL OF PENETRATION DEPTH WHEN DIRECTING HARVESTER TOWARD SCIATIC NOTCH. THE DISTAL ASPECT OF THE SECOND APERTURE IS 5CM FROM THE TIP OF THE TREPHINE.

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[SUGGESTED SIZE- 7MM] 300-070 and 300-07C

- Stay between the inner & outer tables of the ASIS (Anterior Superior Iliac Spinos Process) as you direct the harvester under the rim of the iliac crest.

SUGGESTIONS FOR PROXIMAL TIBIAL METAPHYSEAL HARVESTING

- 1. Make a small (2-3 cm) linear skin incision either just medial (preferred) or lateral to the tibial tubercle.
2. While retracting the skin with a small self-retaining retractor, create a cortical defect.
3. Palpate the joint line and direct the Corex harvester within the metaphysis, remaining 1 to 2 cm below the joint line and closely monitoring the depth of penetration to avoid opposite cortical breach. Bovie coagulate all subcutaneous bleeders and suture the skin.
4. Apply a compressive dressing during the early post-op period.

Manufactured By: Trinity Orthopedics, LLC 8817 Production Avenue, San Diego, CA 92121
Tel: 858-689-4113 fax: 858 689-4115 e-mail: info@trinity-ortho.com
For more information, and instructional videos on product use, go to www.Trinity-Ortho.com

Corex™ II REF 300-090 (0100860413000109) REF 300-09C (0100860413000147) REF 300-070 (0100860413000116) REF 300-07C (0100860413000154)
(Español) Instrucciones de uso

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Fabricado por: Trinity Orthopedics, LLC 8817 Production Avenue, San Diego, CA 92121, EE.UU.